
BOXET 31



Clinical Case Presentation

Inside This Issue: The Encode[®] Complete Restorative System For Fabrication Of Patient Specific Restorations[®]

Case Presentation By: Christopher D. Ramsey, DMD

Volume 1, Issue 2



Treatment Of The Posterior Maxilla With The Encode[®] Complete Restorative System

Clinical Treatment By Christopher D. Ramsey, DMD†

Introduction



Patients desire aesthetic restorations for missing teeth supported by dental implants. With the introduction of CAD/CAM technology for the fabrication of custom abutments, clinicians may be able to provide patients with highly predictable, aesthetic restorations, which support the soft tissue, using a process that mirrors conventional dentistry from ease of use and cost perspectives. Since 2005, it has been possible to fabricate Patient Specific Restorations[®] using Encode

Abutments. Once the surgeon places Encode Healing Abutments at the time of implant placement or at implant uncovering, the restorative dentist need only make a simple impression of the healing abutments. From here, the laboratory technician sends the mounted casts to BIOMET **3***i*.

The newest innovation now includes robotic placement of implant analogs into the master cast in a process called Robocast[™] Technology. This process eliminates the need for an implant level impression, which makes the procedure more efficient and simple, thus saving valuable chairtime for the restorative dentist. Additionally, there is no need to have an inventory of various components. The data from the codes embedded on the occlusal surfaces of Encode Healing Abutments transmits information about the implant connection type, restorative seating surface diameter and hex orientation. The definitive abutments are then designed virtually and the data is transferred to a milling machine for fabrication of definitive Encode Abutments.

The *Clinical Case Presentation* to follow demonstrates a 54-year-old female patient who presented with a failed fixed partial denture supported by teeth Nos. 11 and 13. Tooth No. 13 (maxillary left second bicuspid) was fractured and deemed non-restorable (Figure 1). The patient desired replacement of the failed FPD with fixed restorations. The treatment plan accepted by the patient included extraction of the fractured tooth, implant placement and restoration of teeth Nos. 11, 12, 13 and 14 with individual crowns.

Surgical Treatment

Following administration of local anesthesia, a sulcular incision was made around tooth No. 13, followed by a midcrestal incision in the edentulous area of tooth No. 12. A full thickness mucoperiosteal flap was elevated on the buccal and lingual aspects. The fractured root of tooth No. 13 was carefully extracted using periotomes and the socket was debrided with hand and rotary instruments. A round drill was used to mark the location of the implant in tooth site No. 12. A 2mm diameter twist drill was advanced into the edentulous space (tooth No.12) as well as the extraction site (tooth No. 13) to prepare the osteotomies for implant placement. The drilling protocol continued for each site following the manufacturer's protocol for placement of a 4/5mm x 13mm Full OSSEOTITE XP® (Expanded Platform) Implant for tooth site No. 12 and a 5mm x 11.5mm Full OSSEOTITE® Certain® Implant for tooth site No. 13. Encode Healing Abutments were chosen consistent with the implant platform, appropriate emergence profile and collar height for each implant. The collar heights of the abutments were selected such that the occlusal surfaces of the abutments were at least 1mm supragingival circumferentially. The chosen healing abutments were placed into the internal interfaces of the implants and



BOXET 31











tightened to 20Ncm of torque. The soft tissue flaps were closed around the healing abutments and secured with intermittent resorbable sutures. A periapical radiograph was taken to confirm complete seating of the healing abutments. A prefabricated fixed provisional restoration was placed over the Encode[®] Healing Abutments and was supported by teeth Nos. 11 and 14. The patient was dismissed with post-operative medications and instructions for oral hygiene.

Restorative Treatment

Three months post implant placement, the patient was seen for evaluation. Healing was uneventful. The provisional fixed partial denture was removed from teeth Nos. 11 and 14 revealing healthy soft tissue maturation surrounding the Encode Healing Abutments and tooth preparations (Figure 2). An impression was made of the Encode Healing Abutments and natural tooth preparations. Light body polyvinylsiloxane impression material was syringed around the healing abutments and tooth preparations followed by placement of heavy body impression material in a closed stock tray. The impression tray was seated and the impression material was allowed to set per the manufacturer's instructions. The impression tray was removed and the impression was examined for verification that the entire occlusal surfaces of the Encode Healing Abutments, as well as the soft tissue contours and tooth preparation margins, were recorded (Figure 3). An impression was made of the opposing arch and sent to the commercial laboratory along with a shade selection and occlusal record.

In the laboratory, the impression of the Encode Healing Abutments was poured using a die stone for fabrication of a master cast (Figure 4). The casts were mounted on the recommended articulator (Stratos[™] 100 with Adesso Split Plates, Ivoclar Vivadent, Inc. Amherst, NY) using the occlusal record. The work order was completed and the master casts were sent to BIOMET *3i* for fabrication of definitive Encode Abutments.

After scanning the master cast, the abutments were designed virtually (Figure 5). Using the same data, implant analogs were placed into the cast with a computer driven robotic arm. This process, called Robocast[™] Technology, allows the definitive Encode Abutments to be placed on the master cast (Figure 6) for fabrication of the definitive restorations by the laboratory.

At the insertion appointment, the Encode Healing Abutments were removed from the implants in tooth sites Nos. 12 and 13. The titanium nitride coated definitive Encode Abutments were seated into the internal interfaces of the implants with audible and tactile clicks, which ensured completed seating. The abutments were secured with Gold-Tite® Abutment Screws tightened to 20Ncm of torgue with a torgue driver (Figures 7 & 8). The individual PFM crowns were tried-in, adjusted interproximally and for optimal occlusal contacts in centric and eccentric positions. Verification radiographs were taken to confirm seating of the abutments into the implants and seating of the definitive restorations to the abutments. The definitive restorations for teeth Nos. 11 and 14 were cemented with RelyX™ Vitremer™ Luting Cement (3M ESPE, St. Paul, MN)) and the restorations supported by the implants in tooth sites Nos. 12 and 13 were cemented with Temrex® Cement (Temrex Corporation, Freeport, NY) (Figure 9). The excess cement was carefully removed and a verification radiograph was taken (Figure 10). The patient was dismissed with oral hygiene instructions.

Surgical Colleague: Karina F. Leal, DMD, West Palm Beach, FL Laboratory Colleague: David Haley, ProLab Esthetics, Inc, Naples, FL

†Dr. Ramsey received his dental degree from Temple University School of Dentistry in Philadelphia, Pennsylvania. He has published several articles and is a Clinical Instructor for The Institute for Oral Art and Design in Sarasota, Florida, as well as the Pacific Aesthetic Continuum Live, in Irvine, California. Dr. Ramsey maintains a private practice with a focus on adhesive aesthetic dentistry, in Jupiter, Florida.



Global Headquarters 4555 Riverside Drive Palm Beach Gardens, FL 33410 1-800-342-5454 Outside The U.S.: +1-561-776-6700 Fax: +1-561-776-1272 www.biomet3i.com

SUBSIDIARIES

AUSTRALIA Phone: +61-2-9855-4444 Fax: +61-2-9888-9900

AUSTRIA Phone: +43-(0)6235-200-45 Fax: +43-(0)6235-200-45-9

BELGIUM Phone: +32-2-5410290 Fax: +32-2-5410291

BRAZIL Phone: +55-11-5081-4405 Fax: +55-11-5081-7484

CANADA Phone: +514-956-9843 Fax: +514-956-9844

FRANCE Phone: +33-1-41054343 Fax: +33-1-41054340

GERMANY Phone: +49-721-255177-10 Fax: +49-721-255177-73

IBELAND Phone: +353-1-800-552-752 Fax: +44-1628-820182

MEXICO Phone: +52-55-5679-1619 Fax: +52-55-5684-8098 THE NETHERI ANDS

Phone: +31-(0)78-629-2800 Fax: +31-(0)78-629-2801

NEW ZEALAND Phone: +64-508-122-221 Fax: +64-508-133-331

NORDIC REGION Phone: +46-40-17-6090 Fax: +46-40-17-6099

PORTUGAL Phone: +351-21-000-1647 Fax: +351-21-000-1675

SPAIN Phone: +34-93-470-59-50 Fax: +34-93-372-11-25

SWITZERLAND Phone: +41-44-380-46-46 Fax: +41-44-383-46 55

Phone: +44-800-652-1233

Fax: +44-1628-820182

TURKEY Phone: +90-212-573-6916 Fax: +90-212-662-7635 U.K.

> COLOMBIA *3i* Colombia Phone: +571-612-9362 Fax: +571-620-6412

DISTRIBUTORS ARGENTINA

Dentalmax, SA Phone: +541-1482-71001 Fax: +541-1482-67373

Cybel, SA Phone: +56-2-2321883

Atek Inc. Phone: +86-21-6329-1265 Fax: +86-21-6329-1620

Fax: +56-2-2330176

Sign Up For BIOMET 3i's Electronic Newsletter "BIOMET 3innovations."

Simply Go Online To www.biomet3i.com/signup

CHILE

CHINA

COSTA RICA Implantec S.A. Phone: +506-234-9043 Fax: +506-224-7620

EL SALVADOR Dentimerc SA de CV Phone: +503-263-6350 Fax: +503-263-6676

GREECE Impladend Dental Implants, LLC Phone: +30-2310-501-651 Fax: +30-2310-522-417

H.A. Systems Phone: +972-3-6138777 Fax: +972-3-6138778 HONG KONG Ositek Inc., Ltd. Phone: +852-8121-6601 Fax: +852-3747-3754

ISRAEL

ITALY Biomax, srl. Phone: +39-0444-913410 Fax: +39-0444-913695

JAPAN Implant Innovations Japan Phone: +81-66-868-3012 Fax: +81-66-868-2444

KOREA Jungsan *3i* Corp. Phone: +82-2-516-1808 Fax: +82-2-514-9434 LEBANON *3i* MENA s.a.l. Middle East And North Africa Phone: +961-1-694000

Fax: +961-1-694222 PARAGUAY Andres H. Arce y Cia SRL Phone: +595-21-208185 Fax: +595-21-496291

POI AND Dental Depot Phone: +48-71-341-3091 Fax: +48-71-343-6560

BUSSIA Com-Dental Phone: +7-495-797-6686 Fax: +7-499-242-9567

SINGAPORE Asia Implant Support & Services Phone: +65-6223-2229 Fax: +65-6220-3538

Selective Surgical CC Phone: +27-11-991-7007 Fax: +27-11-672-1391 TAIWAN

Kuo Hwa Dental Suppliers Co., Ltd. Phone: +886-2-2226-1770 Fax: +886-2-2226-8747

THAILAND 3i (Thailand) Co., LTD. Phone: +662-252-6685 Fax: +662-252-6686

UKRAINE Com-Dental Phone: +38-067-7007667 Fax: +38-044-5017117

URUGUAY Pro3implant S.R.L. Phone: +598-2-4034163 Fax: +598-2-4034163



REV A 06/08

SOUTH AFRICA

Certain, Encode, Gold-Tite, OSSEOTITE, OSSEOTITE XP and Patient Specific Restorations are registered trademarks and NanoTite and Corrent, Licolog, out-me, Colo III, Colo III and a later appendix heat on a later and later and a later and later and a later and later and a later an